



TOP GLOVE

TOP QUALITY, TOP EFFICIENCY
GOOD HEALTH, SAFETY FIRST & BE HONEST
THE WORLD'S LARGEST RUBBER GLOVE MANUFACTURER
THE WORLD IS OUR MARKET

SURGICAL GLOVE

POWDERED & POWDER-FREE



Functional Benefits

- Extra strength provides additional protection from surgical debris.
- Fully anatomical design to reduce hand fatigue.
- Softness provides superior comfort and natural fit.
- Micro-roughened surface provides excellent wet and dry grip.
- Easy donning and helps prevent rolling back.
- High strength and elasticity.



surgical glove

SURGICAL GLOVE

POWDERED & POWDER-FREE

Quality Standards

- Conforms to ASTM D3577 and EN455 Standards
- Manufactured under QSR (GMP), ISO9001 : 2008 Quality Management System and ISO 13485:2003
- Using FDA approved absorbable corn starch
- Sterilized by Gamma Ray & ETO
- Bioburden and sterility tested
- Hypoallergenic reducing potential allergic reactions

Glove Sizes

- 5.5, 6.0, 6.5, 7.0, 7.5, 8.0, 8.5, 9.0
- Marked in the check box on the shipping carton

Product Specifications

Type	Gamma ray Sterilized; Powdered or Powder-free
Material	Natural high quality latex
Design & Features	Hand specific; curved fingers; beaded cuff; natural to white, off white to yellow
Storage	The gloves shall maintain their properties when stored in a dry condition at temperature not higher than 30°C
Shelf-life	5 years from date of manufacturing

Physical Dimensions

DIMENSIONS	STANDARDS	
	Top Glove	ASTM D3577
Length (mm)	280 min	245 min (size 5.5) 265 mm (size 6-9)
Width (mm)		
• 5.5	71 ± 5	70 ± 6
• 6.0	77 ± 5	76 ± 6
• 6.5	83 ± 5	83 ± 6
• 7.0	89 ± 5	89 ± 6
• 7.5	95 ± 5	95 ± 6
• 8.0	102 ± 6	102 ± 6
• 8.5	108 ± 6	108 ± 6
• 9.0	115 ± 5	114 ± 6
Thickness - Single Wall (mm)		
• Fingers	0.16 min	0.10 min
• Palm	0.13 min	0.10 min
• Cuff	0.10 min	0.10 min

Physical Properties

Property	ASTM D3577
Elongation at break (%)	
• Before Aging	Min 750
• After Aging	Min 560
Tensile Strength (MPa)	
• Before Aging	Min 24
• After Aging	Min 18



International Quality System Certificate Awarded



GERMANY



U.S.A



AUSTRALIA



Standard Malaysian Glove
Certificate Approval No. 2004



CANADA



QUALITY SYSTEM
MALAYSIA

If you have any enquiries on our products, please contact:

TOP GLOVE GROUP OF COMPANIES
Lot 4969, Jalan Teratai, Batu 6, Off Jalan Meru
41050 Klang, Selangor Darul Ehsan, Malaysia
Tel : +603-3392 1992 / 1905
Fax : +603-3392 1291 / 8410
Email : sales@topglove.com.my
Website : www.topglove.com.my